

# Diskus®



## Tānsi oma ka isi āpahcihtayān oma Diskus®?

1. E-kipahikātek.
2. Ta yotenaman: Waskic anima micihmina peyak ohci kihcihcin ekwa kih mihihcin anima ohci kotak kihcihcin micihmina ita kih micihcin kitakih astek. Mākona kih micihcin esko ta pihtaman ta mahtwitsik.
3. Kih ahtinen anima ka waskawinaman isko ta kih matwitisik.
4. Micimina anima Diskus® otāhcikan awasite kitonik ohci kita pimitaskonaman ekwa ka pakitatāmon.
5. Mewāc e-tahkwātaman anima ka tahkwacikātek, a-iskwātāmo kiyipa **kwayask**. Ekatenā Diskus® otāhcikan. Kāya pakitatāmo mitahtat 10ka matwesihk, ekwa nisihkac pakitatamo.
6. Kipa: kih mihihcinihk ka miciminaman wāyaw ka itinen esko ka kipipayik.

## Ka isi pamihtāyin oma Diskus®

1. Pahkwaha ita kā tahkwātamane-pahkwahk ayānisis ohci ahpo ekinos.
2. Nahastā ita e-pastek, ekā ita e-sāpeyāyak tāsok. Kasihkwew ikamikohk.
3. Diskus® ta kipahikātek oma eka kapatahk; ka yotepitaman piko oma kā wih apacitak.
4. Diskus® ta kih astek ita eka kitahkwatihk, e-kistek ahpo e-cahkāstek ekwa ispimiik e-ayahk kisitew (nawac mena 30°C).
5. Wāpahta anima ita kā wihihkepayik tahtwāw keyapic ka kih apahcihtayin. Kita mihihkopayew ewako neyānan keyapic iskonamanih kotāhcikanihk.