

# Handihaler®



## Tānsi oma ka wi isi apacihtāyan Handihaler®?

Piskina ekwa ka nanistinaman kahkokesiki. Tāwāyi keskipita. Pāskina wewihkinikan ita kā napakāk peyak piko asowacikanis ka nohkwahk. Ispahkenasi anima ka kikamosik ekwa pakitina asowacikans ka pahkitihk. Kanokiskisih: anihi asowacihkanisa. Eka ka kih cahkāsteki ahpo ka sāpopekih.

1. Yotena ka ispāhkihpitaman.
2. Yotena ita kā tahkwātaman e-ispahkepitaman.
3. Peyak asowacikanis posihta it aka kih astek pāmawes āpacihtayin.
4. Kipaha ita kā tahkwātaman kwayask isko ka pihtaman imatwitisik, ekosi e-yotemok kipahikanis.
5. Micimina oma Handihaler® ita kā tahkwātaman e-simatastek ekwa mākon an akā cestahiket ocipāson PEYAKWĀW EKWA PAKITIN.
6. Opime isi pakitatāmo awaste anima ka tahkwātaman.
7. Mahkwāta anima ka tahkwācikātek. Iskwatāmo nisihkēc ekwa itāmihk isko kispānak e-sākaskinecik. Ekatena Handihaler® kitonihk ohci kisik ekā e-yih yihyin isko mitahtat 10 eskwakihtāsowin.
8. Kihtwām tohta 6-7 mena kihtwām ka yih yihyin.
9. Yotena ita kā tahkwātaman ekwa kwatapina anima asowacikanis kā kesapatahk ekwa wepina. Kipaha it a kā tahkwātaman ekwa kipahikan.
10. Kistapāwatā kih cihcihya

# Handihaler®



## Kesi pamitāyin HandiHaler®

Oma HandiHaler ka kih kanāhcihtahk  
pewakwāw peyak pesim ahpo isko e-  
nitawihyihākwāhk.

1. Yotena kipahikan ohpina ita kā tahkwātaman.  
Yoten ka cestahiket ocihpāson ka yotenaman  
wahkitihk.
2. Sepekina HandiHaler® kicih kesowākamik  
nipyi ka kistāpāwek pihko. Ekāweya āpacih  
sepekinikan.
3. PahkwahaHandiHaler® kwayaskekwa e-  
yoteskataman anima kipahikan, it aka  
tahkwātima ekwa wahkitihk ka pahkawayak.