



LUNG SASK

NEBULIZER

Inhaler Delivery Device

HOW TO USE A NEBULIZER

1. Ensure the nebulizer is plugged in & functioning.
2. Attach tubing to compressor air outlet.
3. Unscrew the top of the medication cup.
4. Add the medication to the cup.
5. Screw on the top.
6. Turn on the compressor; ensure the medication is misting.
7. Place either a mask on your face or mouthpiece in your mouth.
8. Breathe in and out until the liquid is all gone from the medication cup.
9. Turn off the nebulizer.



1. Ensure the nebulizer is plugged in & functioning



2. Attach tubing to compressor air outlet



3. Unscrew the top of the medication cup
4. Add the medication to the cup



5. Screw on the top



6. Turn on the compressor; ensure the medication is misting



7. Place either a mask on your face or mouthpiece in your mouth



8. Breathe in and out until the liquid is all gone from the medication cup




9. Turn off the nebulizer

NEBULIZER Inhaler Delivery Device

CARE OF A NEBULIZER

1. Wash the mask and nebulizer medication cup in warm, soapy water.
2. Rinse well and allow them to air-dry before reuse.
3. Replace the tubing and mask/mouthpiece following the brand's recommended timeline.



lungsask.ca



For video instruction:
lungsask.ca/inhaler-resources